BTMF 2/12-C – In vitro diagnosticum

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BTMF 2/12-C – For in vitro diagnostic use

Drugs of abuse in whole blood, lyophilized

Application

This product is intended for use as a control material as part of a laboratory's internal quality assurance process. The control is prepared by spiking human whole blood with drugs of abuse. The concentrations of the drugs and metabolites lie within the expected treatment ranges.

Preparation

Reconstitute the contents by dispensing accurately 2.5 ml of distilled water into the vial. Close the cap and slightly swivel the vial. Mix thoroughly for 20 minutes with a vortex mixer.

To ensure that the solution is homogeneous shortly mix the vial before sampling.

After reconstitution, the sample material is to be handled in the same manner as a patient sample.

Target value

This material was originally prepared for the proficiency test BTMF 2/12 – Drugs of Abuse in whole blood, organized for the GTFCh by ARVECON GmbH. The target values listed are the consensus values obtained from this trial. No restrictions were placed on sample preparation and analytical methodology. Quantitative analyses were performed by the participants using GC/MS and LC/MS(MS). The target values were released by the coordinator of proficiency testing of the GTFCh.

Storage and stability

Storage: + 2° to 8° C

Stability:

- Sealed container: see expiration date given on the package
- After reconstitution, stored tightly capped in the dark: 7 days
- After reconstitution, please extract and assign Cocaine, Benzoyllecgonines and Ecgonine methyl ester on the same day.

Precautions

All materials of human origin should be considered as potentially infectious and treated with the same care as patient specimens.

Each individual blood unit used for the production of the control was tested for the following antigens and antibodies: HBsAG, anti-HIV-1, anti-HIV-2, anti-HBc and anti-HCV and found to be negative.

Lot: 404051295
Order no.: 12WH212C (10 x 2.5 ml)
GTFCH-RE.no: 12-07
Version: 3-201510

Bestandteil / Analyte	Zielwert Target value	Bewertungsgrenzen Control ranges	Konfidenzbereich Confidence range	Einheit Unit
Amphetamin / Amphetamine	88,1	67,7 - 108,5	83,0 - 93,2	ng/ml
Methamphetamin / Methamphetamine	101	77 - 125	97 - 105	ng/ml
MDMA 3,4-Methylenedioxyamphetamine	54	40,4 - 67,6	51,4 - 56,6	ng/ml
MDA 3,4-Methylenedioxyamphetamine	64,2	48,6 - 79,8	60,0 - 68,4	ng/ml
MDE(A) 3,4-Methylenedioxyethylamphetamine	47,4	35,2 - 59,6	45,1 - 49,7	ng/ml
MBDB N-methyl-1-3,4-methylene-dioxyphenyl-2-butamine	59,6	45,0 - 74,2	57,1 - 62,1	ng/ml
Cocain / Cocaine	52,4	39,2 - 65,6	49,7 - 55,1	ng/ml
Benzoyllecgonin / Benzoyllecgonine	137	107 - 167	132 - 142	ng/ml
Ecgoninmethylester / Ecgonine methyl ester	37,4	27,4 - 47,4	34,4 - 40,4	ng/ml
Codein / Codeine	135	105 - 165	129 - 141	ng/ml
Morphin / Morphine	31,7	23,1 - 40,3	29,4 - 34,0	ng/ml
Dihydrocodein / Dihydrocodeine	160	126 - 194	148 - 172	ng/ml
THC Δ^9 -Tetrahydrocannabinol	2,92	1,78 - 4,06	2,58 - 3,26	ng/ml
11-OH-THC 11-OH- Δ^9 -Tetrahydrocannabinol	2,25	1,33 - 3,17	2,06 - 2,44	ng/ml
THC-COOH 11-nor- Δ^9 -THC carbonic acid	43,1	31,9 - 54,3	40,5 - 45,7	ng/ml

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Control ranges

The control ranges were determined by using the standard deviation according to Horwitz. They were determined by the target value and a single standard deviation ($mean \pm SD_{Horwitz}$).

HORWITZ, W. (1982): Evaluation of analytical methods used for regulation of foods and drugs. Anal Chem 54 (1): 67A-76A

Confidence range (of the target value)

The confidence interval indicates the range in which the target value is located with a significance level of 99%.

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