DCT -25% UR – For forensic use only (FFU)

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DCT -25% UR – For forensic use only (FFU)

Drugs of abuse, drugs and alcohol consumption markers in urine, lyophilized

Application

Accuracy control for the internal quality assurance in the laboratory for forensic analytics and for the execution of alcohol and drug screenings. The target values lie +/- 25% above/below the defined cut-off-values of the EWDTs (European Workplace Drug Testing Society). The control material consists of urine on a human basis to which narcotics, drugs, and alcohol consumption markers have been added.

Preparation

Reconstitute the contents by dispensing 2,0 ml of distilled water into the vial, close the cap and mix gently. Allow the vial to sit for 20 minutes before mixing thoroughly with a vortex mixer.

To ensure that the solution is homogeneous shortly mix the vial before sampling.

After reconstitution, the sample material is to be handled in the same manner as a patient sample.

Storage and stability

Storage: + 2° to 8° C

Stability:

- Sealed container: see expiration date given on the package
- After reconstitution, stored tightly capped in the dark: 7 days.
- These molecules degrade in lyophilized state: See the table.

Precautions

All materials of human origin should be considered as potentially infectious and treated with the same care as patient samples.

Lot / Lot:	5401220130
Order no. / Référence:	22UR020A_n2 (10 x 2,0 ml)
Version / Versione:	1 – 202109

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^(*) These molecules degrade in lyophilized state

^(*) According to the annual real-time stability study, these molecules degrade in lyophilized state after 4 years.

^(*) Wrong amount added / Notice: Presence of Haloperidol gravimetric value 650 µg/L – measured 668 µg/L

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Valore atteso $\pm SD_{Horwitz}$

HORWITZ, W. (1982): Evaluation of analytical methods used for regulation of foods and drugs. Anal Chem 54 (1): 67A-76A

Target value

In this lot, the first step taken was to verify the weight of the analytes (target value) using chromatographical methods. Thus, the control ranges relate to the gravimetric values (except for the analytes marked with *³). The control ranges marked with *³ relate to the determined values. The standard deviation was determined using the values for previous lots. The measurements were carried out by laboratories accredited according to DIN ISO/IEC 17025 for forensic toxicology by means of GC/MS and LC/MS(/MS) in line with the guidelines of the GTFCh (Society of Toxicological and Forensic Chemistry in Germany).

The target values will be controlled according to RiliBäk (the guidelines of the German Medical Association) in a control cycle consisting of 15 values. The completed control cycles carried out by the reference laboratories will be made available at www.acq-science.de.

Control range regarding the evaluation criteria of the GTFCh proficiency test

The accepted intervals were determined using the standard deviation according to Horwitz ($SD_{Horwitz}$).

Target value (mean value) $\pm SD_{Horwitz}$

HORWITZ, W. (1982): Evaluation of analytical methods used for regulation of foods and drugs. Anal Chem 54 (1): 67A-76A

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