PPFL SE – In vitro diagnosticum

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Anwendung

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Propofol in serum, lyophilized

Application

This product is intended to be used as a control material as part of laboratory's internal quality assurance processes. The control is prepared by spiking human serum with Propofol. The concentrations of the analytes lie within the expected treatment ranges.

Preparation

Reconstitute the contents by dispensing accurately 1.0 ml of distilled water into the vial. Close the cap and slightly swivel the vial. Mix thoroughly for 20 minutes with a vortex mixer.

To ensure that the solution is homogeneous shortly mix the vial before sampling.

After reconstitution, the sample material is to be handled in the same manner as a patient sample.

Target value

The assigned value was determined by independent laboratories, each accredited to DIN EN 17025. Repeat determinations were carried out daily on 5 days using GC/MS.

Storage and stability

Storage: + 2° to 8° C

Stability:

- Sealed container: see expiration date given on the package
- After reconstitution, stored tightly capped in the dark: 7 days

Precautions

All materials of human origin should be considered as potentially infectious and treated with the same care as patient specimens.

Each individual blood unit used for the production of the control was tested for the following antigens and antibodies: HBsAG, anti-HIV-1, anti-HIV-2, anti-HBc and anti-HCV and found to be negative.

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Control ranges

The control ranges were determined with the help of the standard deviation of the reference laboratory. They result from the target value and the double standard deviation (mean value \pm twice the standard deviation).

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ACQ Science GmbH